

Aufbrennlegierungen		€/g
hochgoldhaltig	<b>GC3 Bio</b> Au 86   Pt 11,7      WAK 14,6	126,90
	<b>GC4 Bio</b> Au 77,0   Pt 19,0   Ag 1,6      WAK 14,5	119,90

Goldgusslegierungen		€/g
hochgoldhaltig	<b>GG 2</b> Au 69,0   Pt 4,0   Pd 6,0   Ag 17,9 16,2      WAK	106,50
	<b>GG 3</b> Au 72,0   Pt 3,3   Ag 13,7   Cu 10,4	107,20
Ag-Pd-Basis	<b>SG 1</b> Pd 39,9   Ag 51,9      WAK 16,6	59,90


Lote		€/g
universal	<b>Lot 1</b> 1050 °C Au 79,0   Pt 3,0   Ag 16,6	109,90
	<b>Lot 9</b> 710 °C Au 72,5   Pt 0,5   Ag 10,0   Cu 3,0	102,50
für NE-Dentallegierungen	<b>Lot 2</b> 1070 °C Au 71,9   Pd 12,0   Ag 5,0   Cu 8,0	106,50
	<b>Lot 5</b> 760 °C Au 73,2   Pt 1,8   Ag 10,6   Cu 3,0	106,50
	<b>Lot 7</b> 930 °C Au 79,8   Ni 15,8   Zn 4,4	109,90

Preisliste Nr. 3 – gültig ab 07.10.2025

# Wichtige Preismitteilung zu Edelmetall-Dentallegierungen



**Combibond**  
Edelmetall-Dentallegierungen



A vertical grey line on the left side of the form area. To its right, there are four horizontal grey lines stacked vertically, serving as a template for text or a list.